Total IgE





Brief Instruction

Preparations:

- 1) Allow microwell plate, reagents and patient sera to come to room temperature
- 2) Dissolve 1 bottle of washing buffer concentrate in 5 liter aqua dest.
- 3) Dissolve 1 tablet p-NPP in one bottle substrate buffer

Assay Procedure:

- 1) Equip corresponding well amount of the microtiter plate according to working protocol. Always fill stripes with empty wells!
- 2) pipet 20 μL calibrators, controls and patient serum, and
 add 200 μL conjugate. Mix well
- 3) Incubate 1 hour at room temperature. Cover plate with lid
- 4) Wash plate 5 times with each 400 μL (automated washer)
- 5) Repeat mixing of substrate solutionpipet 200 µL substrate solution per well
- 6) Incubate 30 minutes at room temperature. Cover plate with lid
- 7) pipet 50 µl stopping solution per well
- 8) Plate measurement at 405 nm (620 nm reference wavelength)
- 9) Data analysis supported by computer software



Bahnhofstrasse 44 48143 Münster, Germany Phone: +49-251-162-53-12

Fax: +49-251-162-53-18

Email: info@bdl-muenster.de

Web: www.bdl-muenster.de